

For Research Use Only

Integrin beta-5 Recombinant antibody, PBS Only (Capture)

Catalog Number: 83514-5-PBS



Basic Information

Catalog Number:

83514-5-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG29783

GenBank Accession Number:

BC006541

GeneID (NCBI):

3693

UNIPROT ID:

P18084

Full Name:

integrin, beta 5

Calculated MW:

88 kDa

Purification Method:

Protein A purification

CloneNo.:

240401H3

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00512-2, Integrin Beta 5 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83514-5-PBS. Detection antibody: 83514-3-PBS. Standard: Ag29783. Range: 0.313-40 ng/mL.