

For Research Use Only

# Mouse PD-1/CD279 Recombinant antibody, PBS Only (Capture)

Catalog Number: 83516-2-PBS



## Basic Information

Catalog Number:

83516-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_008798.2

GeneID (NCBI):

18566

UNIPROT ID:

Q02242

Full Name:

programmed cell death 1

Calculated MW:

32 kDa

Purification Method:

Protein A purification

CloneNo.:

240389F3

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

mouse

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP00496-2, MOUSE PD-1/CD279 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83516-2-PBS. Detection antibody: 83516-5-PBS. Standard: Eg0918. Range: 0.156-20 ng/mL