

For Research Use Only

# CD86 Recombinant antibody, PBS Only (Capture)

Catalog Number: 83523-3-PBS



## Basic Information

Catalog Number:

83523-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_175862

GeneID (NCBI):

942

ENSEMBL Gene ID:

ENSG00000114013

UNIPROT ID:

P42081

Full Name:

CD86 molecule

Calculated MW:

329 aa, 38 kDa

Purification Method:

Protein A purification

CloneNo.:

240429D3

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP00511-1, CD86 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83523-3-PBS. Detection antibody: 83523-1-PBS. Standard: Eg0966. Range: 0.313-40 ng/mL



Cytometric bead array standard curve of MP00511-2, CD86 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83523-3-PBS. Detection antibody: 83523-2-PBS. Standard: Eg0966. Range: 0.313-40 ng/mL