

For Research Use Only

# CEACAM1/CD66a Recombinant antibody, PBS Only (Detector)

Catalog Number: 83524-2-PBS



## Basic Information

Catalog Number:

83524-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC014473

GeneID (NCBI):

634

UNIPROT ID:

P13688

Full Name:

carcinoembryonic antigen-related  
cell adhesion molecule 1 (biliary  
glycoprotein)

Calculated MW:

526 aa, 58 kDa

Purification Method:

Protein A purification

CloneNo.:

240497H8

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP00520-3, CEACAM-1/CD66a Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83524-3-PBS. Detection antibody: 83524-2-PBS. Standard: Eg0314. Range: 0.313-40 ng/mL