

For Research Use Only

CEACAM1/CD66a Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number: 83524-7-PBS



Basic Information

Catalog Number: 83524-7-PBS	GenBank Accession Number: BC014473	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 634	CloneNo.: 253148A12
Isotype: IgG	UNIPROT ID: P13688	
Immunogen Catalog Number: EG7064	Full Name: carcinoembryonic antigen-related cell adhesion molecule 1 (biliary glycoprotein)	
	Calculated MW: 526 aa, 58 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean CEACAM1/CD66a concentration was determined to be 24.1 ng/mL in human serum 1, 10.8 ng/mL in human serum 2.



Cytometric bead array standard curve of MP00520-4, CEACAM1/CD66a Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83524-7-PBS. Detection antibody: 83524-6-PBS. Standard: Eg7064. Range: 0.312-40 ng/mL.