

For Research Use Only

# Rat Adiponectin Recombinant antibody, PBS Only (Capture)

Catalog Number: 83525-1-PBS



## Basic Information

Catalog Number:

83525-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_144744.3

GeneID (NCBI):

246253

UNIPROT ID:

Q8K3R4

Full Name:

adiponectin, C1Q and collagen  
domain containing

Calculated MW:

26 kDa

Purification Method:

Protein A purification

CloneNo.:

240392A7

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

rat

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

# Selected Validation Data



Cytometric bead array standard curve of MP00522-2, RAT Adiponectin Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83525-1-PBS. Detection antibody: 83525-2-PBS. Standard: Eg0872. Range: 0.313-40 ng/mL



Cytometric bead array standard curve of MP00522-3, RAT Adiponectin Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83525-1-PBS. Detection antibody: 83525-3-PBS. Standard: Eg0872. Range: 0.625-80 ng/mL