

For Research Use Only

# LDHC Recombinant antibody, PBS Only (Detector)



Catalog Number: 83526-3-PBS

## Basic Information

Catalog Number:

83526-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_002301

GeneID (NCBI):

3948

UNIPROT ID:

P07864

Full Name:

lactate dehydrogenase C

Calculated MW:

36 kDa

Purification Method:

Protein A purification

CloneNo.:

240502E10

## Applications

Tested Applications:

WB, Indirect ELISA, Cytometric bead array

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



HepG2 cells were subjected to SDS PAGE followed by western blot with 83526-3-RR (LDHC antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 83526-3-PBS in a different storage buffer formulation.



Cytometric bead array standard curve of MP00523-2, LDHC Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83526-1-PBS. Detection antibody: 83526-3-PBS. Standard: SY00609. Range: 39-2500 pg/mL.