

For Research Use Only

# SLC35B4 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83527-2-PBS

## Basic Information

Catalog Number:

83527-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG22669

GenBank Accession Number:

BC008413

GeneID (NCBI):

84912

UNIPROT ID:

Q96950

Full Name:

solute carrier family 35, member B4

Purification Method:

Protein A purification

CloneNo.:

240332G1

## Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP00497-1, SLC35B4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83527-3-PBS. Detection antibody: 83527-2-PBS. Standard: Ag22669. Range: 0.625-80 ng/mL.



Cytometric bead array standard curve of MP00497-2, SLC35B4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83527-1-PBS. Detection antibody: 83527-2-PBS. Standard: Ag22669. Range: 2.5-80 ng/mL.