

For Research Use Only

# CD5 Recombinant antibody, PBS Only (Capture/Detector)



Catalog Number: 83543-5-PBS

## Basic Information

Catalog Number:

83543-5-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_014207.4

GeneID (NCBI):

921

ENSEMBL Gene ID:

ENSG00000110448

UNIPROT ID:

P06127

Full Name:

CD5 molecule

Calculated MW:

55 kDa

Purification Method:

Protein A purification

CloneNo.:

240428E10

## Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP00525-1, CD5 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83543-5-PBS. Detection antibody: 83543-4-PBS. Standard: Eg1119. Range: 0.313-40 ng/mL



Cytometric bead array standard curve of MP00525-3, CD5 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83543-2-PBS. Detection antibody: 83543-5-PBS. Standard: Eg1119. Range: 0.313-40 ng/mL