

For Research Use Only

# Mouse IL-2 Recombinant antibody, PBS Only (Capture/Detector)



Catalog Number: 83547-2-PBS

## Basic Information

Catalog Number:

83547-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_008366

GeneID (NCBI):

16183

UNIPROT ID:

P04351

Full Name:

interleukin 2

Calculated MW:

19 kDa

Purification Method:

Protein A purification

CloneNo.:

240358D6

## Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Mouse

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP00533-1, MOUSE IL-2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83547-1-PBS. Detection antibody: 83547-2-PBS. Standard: Eg0762. Range: 0.313-40 ng/mL



Cytometric bead array standard curve of MP00533-2, MOUSE IL-2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83547-2-PBS. Detection antibody: 83547-3-PBS. Standard: Eg0762. Range: 0.313-40 ng/mL