

For Research Use Only

Mouse Thrombopoietin Recombinant antibody, PBS Only (Capture)

Catalog Number: 83549-1-PBS



Basic Information

Catalog Number:

83549-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG0678

GenBank Accession Number:

NM_001173505.1

GeneID (NCBI):

21832

UNIPROT ID:

P40226

Full Name:

thrombopoietin

Calculated MW:

38KD

Purification Method:

Protein A purification

CloneNo.:

240203A2

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00544-1, MOUSE Thrombopoietin Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83549-1-PBS. Detection antibody: 83549-2-PBS. Standard: Eg0678. Range: 0.156-20 ng/mL



Cytometric bead array standard curve of MP00544-3, MOUSE Thrombopoietin Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83549-1-PBS. Detection antibody: 83549-4-PBS. Standard: Eg0678. Range: 0.156-20 ng/mL