

For Research Use Only

Mouse Chemerin Recombinant antibody, PBS Only (Detector)



Catalog Number: 83550-3-PBS

Basic Information

Catalog Number:

83550-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_001347168.1

GeneID (NCBI):

71660

UNIPROT ID:

Q9DD06

Full Name:

retinoic acid receptor responder
(tazarotene induced) 2

Calculated MW:

18 kDa

Purification Method:

Protein A purification

CloneNo.:

240355B12

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00537-2, MOUSE Chemerin Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83550-2-PBS. Detection antibody: 83550-3-PBS. Standard: Eg0875. Range: 0.313-40 ng/mL