

For Research Use Only

CD138/Syndecan-1 Recombinant antibody, PBS Only (Detector)

Catalog Number: 83551-4-PBS



Basic Information

Catalog Number:

83551-4-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC008765

GeneID (NCBI):

6382

UNIPROT ID:

P18827

Full Name:

syndecan 1

Calculated MW:

32.5 kDa

Purification Method:

Protein A purification

CloneNo.:

240422C12

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00527-1, CD138/Syndecan-1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83551-2-PBS. Detection antibody: 83551-4-PBS. Standard: Eg0132. Range: 0.313-40 ng/mL



Cytometric bead array standard curve of MP00527-2, CD138/Syndecan-1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83551-3-PBS. Detection antibody: 83551-4-PBS. Standard: Eg0132. Range: 0.313-40 ng/mL