

For Research Use Only

TMEM53 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83554-1-PBS

Basic Information

Catalog Number: 83554-1-PBS	GenBank Accession Number: BC064520	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 79639	CloneNo.: 240511E2
Source: Rabbit	UNIPROT ID: Q6P2H8	
Isotype: IgG	Full Name: transmembrane protein 53	
Immunogen Catalog Number: AG21321	Calculated MW: 277 aa, 32 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00539-1, TMEM53 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83554-2-PBS. Detection antibody: 83554-1-PBS. Standard: Ag21321. Range: Ag21321.