

For Research Use Only

NLRP6 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83555-4-PBS

Basic Information

Catalog Number: 83555-4-PBS	GenBank Accession Number: NM_138329.2	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 171389	CloneNo.: 240543H4
Source: Rabbit	UNIPROT ID: P59044	
Isotype: IgG	Full Name: NLR family, pyrin domain containing 6	
Immunogen Catalog Number: AG34516	Calculated MW: 99 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Cytometric bead array standard curve of MP00540-1, NLRP6 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83555-4-PBS. Detection antibody: 83555-1-PBS. Standard: Ag34516. Range: 0.78-100 ng/mL



Cytometric bead array standard curve of MP00540-2, NLRP6 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83555-4-PBS. Detection antibody: 83555-3-PBS. Standard: Ag34516. Range: 0.78-100 ng/mL



Cytometric bead array standard curve of MP00540-3, NLRP6 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83555-4-PBS. Detection antibody: 83555-2-PBS. Standard: Ag34516. Range: 0.78-100 ng/mL