

For Research Use Only

# MYO1G Recombinant antibody, PBS Only (Capture)

Catalog Number: 83556-2-PBS



## Basic Information

Catalog Number:

83556-2-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG22698

GenBank Accession Number:

BC113544

GeneID (NCBI):

64005

UNIPROT ID:

B0I1T2

Full Name:

myosin IG

Calculated MW:

1018 aa, 116 kDa

Purification Method:

Protein A purification

CloneNo.:

240503C7

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP00541-2, MYO1G Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83556-2-PBS. Detection antibody: 83556-1-PBS. Standard: Ag22698. Range: 0.78-100 ng/mL