

For Research Use Only

# GPR161 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83557-1-PBS

## Basic Information

Catalog Number:	83557-1-PBS	GenBank Accession Number:	BC028163	Purification Method:	Protein A purification
Size:	1 mg/ml	GeneID (NCBI):	23432	CloneNo.:	240522A8
Source:	Rabbit	UNIPROT ID:	Q8N6U8		
Isotype:	IgG	Full Name:	G protein-coupled receptor 161		
Immunogen Catalog Number:	AG30998	Calculated MW:	529 aa, 59 kDa		

## Applications

**Tested Applications:**  
Indirect ELISA, Cytometric bead array

**Species Specificity:**  
Human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP00542-3, GPR161 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83557-4-PBS. Detection antibody: 83557-1-PBS. Standard: Ag30998. Range: 0.313- 20 ng/mL