

For Research Use Only

# GPR161 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83557-4-PBS

## Basic Information

Catalog Number: 83557-4-PBS	GenBank Accession Number: BC028163	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 23432	CloneNo.: 240522D8
Source: Rabbit	UNIPROT ID: Q8N6U8	
Isotype: IgG	Full Name: G protein-coupled receptor 161	
Immunogen Catalog Number: AG30998	Calculated MW: 529 aa, 59 kDa	

## Applications

Tested Applications:  
Indirect ELISA, Cytometric bead array

Species Specificity:  
Human

## Background Information

### Storage

Storage:  
Store at -80°C.  
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP00542-2, GPR161 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83557-4-PBS. Detection antibody: 83557-3-PBS. Standard: Ag30998. Range: 0.156- 20 ng/mL



Cytometric bead array standard curve of MP00542-3, GPR161 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83557-4-PBS. Detection antibody: 83557-1-PBS. Standard: Ag30998. Range: 0.313- 20 ng/mL