

For Research Use Only

S100A14 Recombinant antibody, PBS Only (Capture/Detector)



Catalog Number: 83559-4-PBS

Basic Information

Catalog Number:	83559-4-PBS	GenBank Accession Number:	BC005019	Purification Method:	Protein A purification
Size:	1 mg/ml	GeneID (NCBI):	57402	CloneNo.:	240528G11
Source:	Rabbit	UNIPROT ID:	Q9HCY8		
Isotype:	IgG	Full Name:	S100 calcium binding protein A14		
Immunogen Catalog Number:	AG0757	Calculated MW:	12 kDa		

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00547-2, S100A14 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83559-1-PBS. Detection antibody: 83559-4-PBS. Standard: Ag0757. Range: 0.781-100 ng/mL



Cytometric bead array standard curve of MP00547-3, S100A14 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83559-3-PBS. Detection antibody: 83559-4-PBS. Standard: Ag0757. Range: 0.781-100 ng/mL