

For Research Use Only

ARID3B Recombinant antibody, PBS Only (Capture/Detector)



Catalog Number: 83560-2-PBS

Basic Information

Catalog Number:

83560-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG10545

GenBank Accession Number:

BC041792

GeneID (NCBI):

10620

UNIPROT ID:

Q8IVW6

Full Name:

AT rich interactive domain 3B (BRIGHT-like)

Calculated MW:

61 kDa

Purification Method:

Protein A purification

CloneNo.:

240466F5

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00548-2, ARID3B Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83560-2-PBS. Detection antibody: 83560-3-PBS. Standard: Ag10545. Range: 0.781-100 ng/mL.



Cytometric bead array standard curve of MP00548-3, ARID3B Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83560-1-PBS. Detection antibody: 83560-2-PBS. Standard: Ag10545. Range: 0.781-100 ng/mL.