

For Research Use Only

SYTL4 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83564-1-PBS

Basic Information

Catalog Number: 83564-1-PBS	GenBank Accession Number: BC014913	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 94121	CloneNo.: 240556C2
Source: Rabbit	UNIPROT ID: Q96C24	
Isotype: IgG	Full Name: synaptotagmin-like 4	
Immunogen Catalog Number: AG2772	Calculated MW: 671 aa, 76 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00538-2, SYTL4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83564-2-PBS. Detection antibody: 83564-1-PBS. Standard: Ag2772. Range: 0.625-80 ng/mL



Cytometric bead array standard curve of MP00538-3, SYTL4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83564-4-PBS. Detection antibody: 83564-1-PBS. Standard: Ag2772. Range: 1.25-80 ng/mL