

For Research Use Only

SYTL4 Recombinant antibody, PBS Only (Capture/Detector)

Catalog Number: 83564-4-PBS



Basic Information

Catalog Number:

83564-4-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2772

GenBank Accession Number:

BC014913

GeneID (NCBI):

94121

UNIPROT ID:

Q96C24

Full Name:

synaptotagmin-like 4

Calculated MW:

671 aa, 76 kDa

Purification Method:

Protein A purification

CloneNo.:

240556F9

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00538-1, SYTL4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83564-3-PBS. Detection antibody: 83564-4-PBS. Standard: Ag2772. Range: 0.625-80 ng/mL



Cytometric bead array standard curve of MP00538-3, SYTL4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83564-4-PBS. Detection antibody: 83564-1-PBS. Standard: Ag2772. Range: 1.25-80 ng/mL