

For Research Use Only

SYTL4 Recombinant antibody, PBS Only



Catalog Number: 83564-5-PBS

Basic Information

Catalog Number:

83564-5-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2772

GenBank Accession Number:

BC014913

GeneID (NCBI):

94121

UNIPROT ID:

Q96C24

Full Name:

synaptotagmin-like 4

Calculated MW:

671 aa, 76 kDa

Purification Method:

Protein A purification

CloneNo.:

240556E2

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 83564-5-RR (SYTL4 antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 83564-5-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 83564-5-RR against Human SYTL4 were performed. The affinity constant is below 1 pM.