

For Research Use Only

TMEM16A/DOG1 Recombinant antibody, PBS Only (Detector)

Catalog Number: 83567-4-PBS



Basic Information

Catalog Number:

83567-4-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3320

GenBank Accession Number:

BC027590

GeneID (NCBI):

55107

UNIPROT ID:

Q5XXA6

Full Name:

anoctamin 1, calcium activated chloride channel

Calculated MW:

986 aa, 114 kDa

Purification Method:

Protein A purification

CloneNo.:

241464G10

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01115-2, ANO1/TMEM16A/DOG1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83567-2-PBS. Detection antibody: 83567-4-PBS. Standard: Ag3320. Range: 0.313-20 ng/mL



Biolayer interferometry (BLI) kinetic assay of 83567-4-PBS against Human TMEM16A/DOG1 was performed. The affinity constant is below 1 pM.