

For Research Use Only

# DMC1 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83568-1-PBS

## Basic Information

|                                     |   |  |
|-------------------------------------|---|--|
| Catalog Number:<br>83568-1-PBS      | GenBank Accession Number:<br>BC035658   | Purification Method:<br>Protein A purification |
| Size:<br>1 mg/ml                    | GeneID (NCBI):<br>11144   | CloneNo.:<br>240463A7                          |
| Source:<br>Rabbit                   | UNIPROT ID:<br>Q14565   |  |
| Isotype:<br>IgG                     | Full Name:<br>DMC1 dosage suppressor of mck1<br>homolog, meiosis-specific<br>homologous recombination (yeast) |  |
| Immunogen Catalog Number:<br>AG4653 | Calculated MW:<br>165 aa, 18 kDa, 38 kDa  |  |

## Applications

Tested Applications:  
Indirect ELISA, Cytometric bead array

Species Specificity:  
Human

## Background Information

### Storage

Storage:  
Store at -80°C.  
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP00558-3, DMC1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83568-5-PBS. Detection antibody: 83568-1-PBS. Standard: Ag4653. Range: 0.625-80 ng/mL