

For Research Use Only

DMC1 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83568-4-PBS

Basic Information

Catalog Number: 83568-4-PBS	GenBank Accession Number: BC035658	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 11144	CloneNo.: 240463B11
Source: Rabbit	UNIPROT ID: Q14565	
Isotype: IgG	Full Name: DMC1 dosage suppressor of mck1 homolog, meiosis-specific homologous recombination (yeast)	
Immunogen Catalog Number: AG4653	Calculated MW: 165 aa, 18 kDa, 38 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00558-1, DMC1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83568-2-PBS. Detection antibody: 83568-4-PBS. Standard: Ag4653. Range: 0.625-80 ng/mL



Cytometric bead array standard curve of MP00558-2, DMC1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83568-3-PBS. Detection antibody: 83568-4-PBS. Standard: Ag4653. Range: 0.625-80 ng/mL