

For Research Use Only

# ANAPC10 Recombinant antibody, PBS Only (Capture/Detector)



Catalog Number: 83570-2-PBS

## Basic Information

Catalog Number:

83570-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG8487

GenBank Accession Number:

BC005217

GeneID (NCBI):

10393

UNIPROT ID:

Q9UM13

Full Name:

anaphase promoting complex subunit

10

Calculated MW:

185 aa, 21 kDa

Purification Method:

Protein A purification

CloneNo.:

240492F12

## Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP00563-1, ANAPC10 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83570-2-PBS. Detection antibody: 83570-3-PBS. Standard: Ag8487. Range: 0.313-40 ng/mL



Cytometric bead array standard curve of MP00563-2, ANAPC10 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83570-1-PBS. Detection antibody: 83570-2-PBS. Standard: Ag8487. Range: 0.313-40 ng/mL