

For Research Use Only

Mouse VEGFA164 Recombinant antibody, PBS Only (Detector)

Catalog Number: 83572-2-PBS



Basic Information

Catalog Number:

83572-2-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG0277

GenBank Accession Number:

NM_001287057

GeneID (NCBI):

22339

UNIPROT ID:

Q00731-2

Full Name:

vascular endothelial growth factor A

Calculated MW:

22kd

Purification Method:

Protein A purification

CloneNo.:

240415G2

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00561-1, MOUSE VEGF Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83572-4-PBS. Detection antibody: 83572-2-PBS. Standard: Eg0277. Range: 0.625-80 ng/mL



Cytometric bead array standard curve of MP00561-2, MOUSE VEGF Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83572-5-PBS. Detection antibody: 83572-2-PBS. Standard: Eg0277. Range: 0.625-80 ng/mL