

For Research Use Only

CD320 Recombinant antibody, PBS Only (Capture/Detector)



Catalog Number: 83574-2-PBS

Basic Information

Catalog Number:

83574-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG21445

GenBank Accession Number:

BC000668

GeneID (NCBI):

51293

UNIPROT ID:

Q9NPF0

Full Name:

CD320 molecule

Calculated MW:

29 kDa

Purification Method:

Protein A purification

CloneNo.:

240500H5

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00564-2, CD320 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83574-2-PBS. Detection antibody: 83574-4-PBS. Standard: Ag21445. Range: 0.625-40 ng/mL



Cytometric bead array standard curve of MP00564-3, CD320 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83574-3-PBS. Detection antibody: 83574-2-PBS. Standard: Ag21445. Range: 1.25-40 ng/mL