

For Research Use Only



HIP14; ZDHHC17 Recombinant antibody, PBS Only (Capture)

Catalog Number: 83575-1-PBS

Basic Information

Catalog Number:	83575-1-PBS	GenBank Accession Number:	BC050324	Purification Method:	Protein A purification
Size:	1 mg/ml	GeneID (NCBI):	23390	CloneNo.:	240515A3
Source:	Rabbit	UNIPROT ID:	Q8IUH5		
Isotype:	IgG	Full Name:	zinc finger, DHHC-type containing 17		
Immunogen Catalog Number:	AG32654	Calculated MW:	73 kDa		

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00566-2, HIP14; ZDHHC17 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83575-1-PBS. Detection antibody: 83575-3-PBS. Standard: Ag32654. Range: 0.313-40 ng/mL



Cytometric bead array standard curve of MP00566-3, HIP14; ZDHHC17 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83575-1-PBS. Detection antibody: 83575-2-PBS. Standard: Ag32654. Range: 0.625-80 ng/mL