

For Research Use Only

HIP14; ZDHHC17 Recombinant antibody, PBS Only (Capture/Detector)

Catalog Number: **83575-2-PBS**

Basic Information

Catalog Number: 83575-2-PBS	GenBank Accession Number: BC050324	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 23390	CloneNo.: 240515C1
Source: Rabbit	UNIPROT ID: Q8IUH5	
Isotype: IgG	Full Name: zinc finger, DHHC-type containing 17	
Immunogen Catalog Number: AG32654	Calculated MW: 73 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00566-1, HIP14; ZDHHC17 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83575-2-PBS. Detection antibody: 83575-4-PBS. Standard: Ag32654. Range: 0.313-40 ng/mL.



Cytometric bead array standard curve of MP00566-3, HIP14; ZDHHC17 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83575-1-PBS. Detection antibody: 83575-2-PBS. Standard: Ag32654. Range: 0.625-80 ng/mL.