

For Research Use Only

HBQ1 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83576-1-PBS

Basic Information

Catalog Number:	83576-1-PBS	GenBank Accession Number:	BC056686	Purification Method:	Protein A purification
Size:	1 mg/ml	GeneID (NCBI):	3049	CloneNo.:	240527A10
Source:	Rabbit	UNIPROT ID:	P09105		
Isotype:	IgG	Full Name:	hemoglobin, theta 1		
Immunogen Catalog Number:	AG10105	Calculated MW:	142 aa, 16 kDa		

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00567-1, HBQ1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83576-1-PBS. Detection antibody: 83576-2-PBS. Standard: Ag10105. Range: 0.781-100 ng/mL.