

For Research Use Only

HBQ1 Recombinant antibody, PBS Only (Capture/Detector)



Catalog Number: 83576-2-PBS

Basic Information

Catalog Number:

83576-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG10105

GenBank Accession Number:

BC056686

GeneID (NCBI):

3049

UNIPROT ID:

P09105

Full Name:

hemoglobin, theta 1

Calculated MW:

142 aa, 16 kDa

Purification Method:

Protein A purification

CloneNo.:

240527G12

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00567-1, HBQ1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83576-1-PBS. Detection antibody: 83576-2-PBS. Standard: Ag10105. Range: 0.781-100 ng/mL.



Cytometric bead array standard curve of MP00567-2, HBQ1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83576-2-PBS. Detection antibody: 83576-3-PBS. Standard: Ag10105. Range: 0.781-100 ng/mL.