

For Research Use Only

# HBQ1 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83576-3-PBS

## Basic Information

Catalog Number:

83576-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG10105

GenBank Accession Number:

BC056686

GeneID (NCBI):

3049

UNIPROT ID:

P09105

Full Name:

hemoglobin, theta 1

Calculated MW:

142 aa, 16 kDa

Purification Method:

Protein A purification

CloneNo.:

240527C2

## Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP00567-2, HBQ1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83576-2-PBS. Detection antibody: 83576-3-PBS. Standard: Ag10105. Range: 0.781-100 ng/mL.