

For Research Use Only

GAR1 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83577-3-PBS

Basic Information

Catalog Number:	83577-3-PBS	GenBank Accession Number:	BC003413	Purification Method:	Protein A purification
Size:	1 mg/ml	GeneID (NCBI):	54433	CloneNo.:	240537G10
Source:	Rabbit	UNIPROT ID:	Q9NY12		
Isotype:	IgG	Full Name:	GAR1 ribonucleoprotein homolog (yeast)		
Immunogen Catalog Number:	AG2282	Calculated MW:	217 aa, 22 kDa		

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00568-2, GAR1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83577-3-PBS. Detection antibody: 83577-1-PBS. Standard: Ag2282. Range: 0.781-100 ng/mL.