

For Research Use Only

TRIM15 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83582-2-PBS

Basic Information

Catalog Number: 83582-2-PBS	GenBank Accession Number: BC038585	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 89870	CloneNo.: 240565D4
Source: Rabbit	UNIPROT ID: Q9C019	
Isotype: IgG	Full Name: tripartite motif-containing 15	
Immunogen Catalog Number: AG4533	Calculated MW: 465 aa, 52 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00555-2, TRIM15 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83582-2-PBS. Detection antibody: 83582-3-PBS. Standard: Ag4533. Range: 0.78-100 ng/mL



Cytometric bead array standard curve of MP00555-3, TRIM15 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83582-2-PBS. Detection antibody: 83582-4-PBS. Standard: Ag4533. Range: 0.78-100 ng/mL