

For Research Use Only

MCM9 Recombinant antibody, PBS Only



Catalog Number: 83587-1-PBS

Basic Information

Catalog Number:

83587-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG4108

GenBank Accession Number:

BC031658

GeneID (NCBI):

254394

UNIPROT ID:

Q9NXL9

Full Name:

minichromosome maintenance
complex component 9

Calculated MW:

1143 aa, 127 kDa

Purification Method:

Protein A purification

CloneNo.:

240252E9

Applications

Tested Applications:

WB, ELISA

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 83587-1-RR (MCM9 antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 83587-1-PBS in a different storage buffer formulation.