

For Research Use Only

Drebrin Recombinant antibody, PBS Only (Capture)



Catalog Number: 83591-3-PBS

Basic Information

Catalog Number: 83591-3-PBS	GenBank Accession Number: BC000283	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 1627	CloneNo.: 240599E10
Source: Rabbit	ENSEMBL Gene ID: ENSG00000113758	
Isotype: IgG	UNIPROT ID: Q16643	
Immunogen Catalog Number: AG0371	Full Name: drebrin 1	
	Calculated MW: 71 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00585-1, Drebler Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83591-3-PBS. Detection antibody: 83591-4-PBS. Standard: Ag0371. Range: 0.78-100 ng/mL