

For Research Use Only

Neurogranin Recombinant antibody, PBS Only (Detector)



Catalog Number: 83593-2-PBS

Basic Information

Catalog Number: 83593-2-PBS	GenBank Accession Number: BC002835	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 4900	CloneNo.: 240526F12
Source: Rabbit	UNIPROT ID: Q92686	
Isotype: IgG	Full Name: neurogranin (protein kinase C substrate, RC3)	
Immunogen Catalog Number: AG0676	Calculated MW: 7.6 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00606-2, NRG1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83593-3-PBS. Detection antibody: 83593-2-PBS. Standard: Ag0676. Range: 0.156-20 ng/mL