

For Research Use Only

KMT2C Recombinant antibody, PBS Only (Capture)

Catalog Number: 83595-3-PBS



Basic Information

Catalog Number:

83595-3-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG28875

GenBank Accession Number:

NM_170606

GeneID (NCBI):

58508

UNIPROT ID:

Q8NEZ4

Full Name:

myeloid/lymphoid or mixed-lineage leukemia 3

Calculated MW:

541 kDa

Purification Method:

Protein A purification

CloneNo.:

240542E4

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00610-2, KMT2C Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83595-3-PBS. Detection antibody: 83595-2-PBS. Standard: Ag28875. Range: 0.313-40 ng/mL



Cytometric bead array standard curve of MP00610-3, KMT2C Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83595-3-PBS. Detection antibody: 83595-1-PBS. Standard: Ag28875. Range: 0.313-40 ng/mL