

For Research Use Only

BARHL2 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83596-1-PBS

Basic Information

| | | |
|--------------------------------------|---------------------------------------|--|
| Catalog Number: 83596-1-PBS | GenBank Accession Number: BC126439 | Purification Method: Protein A purification |
| Size: 1 mg/ml | GeneID (NCBI): 343472 | CloneNo.: 240558B5 |
| Source: Rabbit | UNIPROT ID: Q9NY43 | |
| Isotype: IgG | Full Name: BarH-like homeobox 2 | |
| Immunogen Catalog Number: AG21156 | Calculated MW: 387 aa, 42 kDa | |

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00582-1, BARHL2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83596-1-PBS. Detection antibody: 83596-2-PBS. Standard: Ag21156. Range: 0.313-40 ng/mL



Cytometric bead array standard curve of MP00582-2, BARHL2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83596-1-PBS. Detection antibody: 83596-4-PBS. Standard: Ag21156. Range: 0.313-40 ng/mL