

For Research Use Only

# SLC39A14/ZIP-14 Recombinant antibody, PBS Only (Capture)

Catalog Number: 83597-1-PBS



## Basic Information

Catalog Number:

83597-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG24210

GenBank Accession Number:

BC015770

GeneID (NCBI):

23516

UNIPROT ID:

Q15043

Full Name:

solute carrier family 39 (zinc transporter), member 14

Calculated MW:

54 kDa

Purification Method:

Protein A purification

CloneNo.:

240577A6

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

# Selected Validation Data



Cytometric bead array standard curve of MP00583-2, SLC39A14/ZIP-14 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83597-1-PBS. Detection antibody: 83597-2-PBS. Standard: Ag24210. Range: 0.78-100 ng/mL.