

For Research Use Only

# Mouse Keap1 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83599-2-PBS

## Basic Information

Catalog Number: 83599-2-PBS	GenBank Accession Number: NM_001110305	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 50868	CloneNo.: 240549G8
Source: Rabbit	UNIPROT ID: Q9Z2X8	
Isotype: IgG	Full Name: kelch-like ECH-associated protein 1	
Immunogen Catalog Number: AG35035	Calculated MW: 70 kDa	

## Applications

Tested Applications:  
Indirect ELISA, Cytometric bead array

Species Specificity:  
Mouse

## Background Information

### Storage

Storage:  
Store at -80°C.  
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP00586-2, MOUSE Keap1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83599-2-PBS. Detection antibody: 83599-3-PBS. Standard: Ag35035. Range: 0.625-80 ng/mL.



Cytometric bead array standard curve of MP00586-3, MOUSE Keap1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83599-2-PBS. Detection antibody: 83599-4-PBS. Standard: Ag35035. Range: 0.625-80 ng/mL.