

For Research Use Only

# PD-L1/CD274 Recombinant antibody, PBS Only (Detector)

Catalog Number: 83600-3-PBS



## Basic Information

Catalog Number:

83600-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC074984

GeneID (NCBI):

29126

UNIPROT ID:

Q9NZQ7

Full Name:

CD274 molecule

Calculated MW:

290 aa, 33 kDa

Purification Method:

Protein A purification

CloneNo.:

240721C3

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP00601-2, PD-L1/CD274 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83600-1-PBS. Detection antibody: 83600-3-PBS. Standard: Eg1114. Range: 0.156-20 ng/mL.