

For Research Use Only

TMPRSS4 Recombinant antibody, PBS Only (Capture)

Catalog Number: 83604-1-PBS



Basic Information

Catalog Number: 83604-1-PBS	GenBank Accession Number: BC011703	Purification Method: Protein A purification
Concentration: 1 mg/ml	GeneID (NCBI): 56649	CloneNo.: 230282E6
Source: Rabbit	UNIPROT ID: Q9NRS4	
Isotype: IgG	Full Name: transmembrane protease, serine 4	
Immunogen Catalog Number: AG1804	Calculated MW: 48 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00581-1, TMPRSS4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83604-1-PBS. Detection antibody: 83604-3-PBS. Standard: Ag1804. Range: 0.391-50 ng/mL



Cytometric bead array standard curve of MP00581-2, TMPRSS4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83604-1-PBS. Detection antibody: 83604-2-PBS. Standard: Ag1804. Range: 0.781-100 ng/mL