

For Research Use Only

# PIK3R2 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83606-3-PBS

## Basic Information

Catalog Number:

83606-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG14367

GenBank Accession Number:

BC014170

GeneID (NCBI):

5296

UNIPROT ID:

O00459

Full Name:

phosphoinositide-3-kinase,  
regulatory subunit 2 (beta)

Calculated MW:

728 aa, 82 kDa

Purification Method:

Protein A purification

CloneNo.:

240531D11

## Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP00587-3, PIK3R2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83606-2-PBS. Detection antibody: 83606-3-PBS. Standard: Ag14367. Range: 0.625-80 ng/mL