

For Research Use Only

ATM Recombinant antibody, PBS Only (Capture)

Catalog Number: 83608-1-PBS



Basic Information

Catalog Number:

83608-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG25920

GenBank Accession Number:

BC137169

GeneID (NCBI):

472

UNIPROT ID:

Q13315

Full Name:

ataxia telangiectasia mutated

Purification Method:

Protein A purification

CloneNo.:

240560B1

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00592-1, ATM Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83608-1-PBS. Detection antibody: 83608-4-PBS. Standard: Ag25920. Range: 0.781-100 ng/mL.



Cytometric bead array standard curve of MP00592-2, ATM Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83608-1-PBS. Detection antibody: 83608-3-PBS. Standard: Ag25920. Range: 0.781-100 ng/mL.