

For Research Use Only

Rat VEGFA Recombinant antibody, PBS Only

Catalog Number: 83610-3-PBS



Basic Information

Catalog Number:

83610-3-PBS

Concentration:

1000 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG0941

GenBank Accession Number:

NM_001287108.1

GeneID (NCBI):

83785

UNIPROT ID:

P16612-2

Full Name:

vascular endothelial growth factor A

Calculated MW:

25 kDa

Purification Method:

Protein A purification

CloneNo.:

240573H8

Applications

Tested Applications:

Indirect ELISA

Species Specificity:

rat

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

