

For Research Use Only

AP1M1 Recombinant antibody, PBS Only (Capture)

Catalog Number: 83613-2-PBS



Basic Information

Catalog Number:

83613-2-PBS

Concentration:

1000 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2757

GenBank Accession Number:

BC017469

GeneID (NCBI):

8907

UNIPROT ID:

Q9BXS5

Full Name:

adaptor-related protein complex 1,
mu 1 subunit

Calculated MW:

423 aa, 49 kDa

Purification Method:

Protein A purification

CloneNo.:

240553A9

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00588-1, AP1M1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83613-2-PBS. Detection antibody: 83613-5-PBS. Standard: Ag2757. Range: 0.781-100 ng/mL.



Cytometric bead array standard curve of MP00588-2, AP1M1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83613-2-PBS. Detection antibody: 83613-4-PBS. Standard: Ag2757. Range: 3.125-100 ng/mL.