

For Research Use Only

# AP1M1 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83613-5-PBS

## Basic Information

Catalog Number: 83613-5-PBS	GenBank Accession Number: BC017469	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 8907	CloneNo.: 240553G9
Source: Rabbit	UNIPROT ID: Q9BX55	
Isotype: IgG	Full Name: adaptor-related protein complex 1, mu 1 subunit	
Immunogen Catalog Number: AG2757	Calculated MW: 423 aa, 49 kDa	

## Applications

Tested Applications:  
Indirect ELISA, Cytometric bead array

Species Specificity:  
Human

## Background Information

### Storage

Storage:  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

# Selected Validation Data



Cytometric bead array standard curve of MP00588-1, AP1M1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83613-2-PBS. Detection antibody: 83613-5-PBS. Standard: Ag2757. Range: 0.781-100 ng/mL.